# visomat

## comfort XXL

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## Table of contents

<b>A</b>	<b>Intended use</b>	34
<b>B</b>	<b>Safety instructions</b>	
	1. Important patient information	35
	2. Important technical information	36
<b>C</b>	<b>Operating the device</b>	
	1. Device description	39
	2. Display	40
	3. Important usage information	41
	4. Initial operation of the device	42
	5. Inserting/replacing batteries, battery icon	42
	6. Activating and setting the date/time	43
	7. Attaching the cuff	44
	8. Body position during the measurement	45
	9. Selecting the user	46
	10. Measuring blood pressure	46
	11. Classification of the measured values (WHO)	47
	12. Irregular pulses and cardiac arrhythmia	48
	13. Using the memory	48
<b>D</b>	<b>What you should know about blood pressure</b>	
	1. The systolic and diastolic blood pressure value	50
	2. Why you measure different values	50
	3. Why should you measure blood pressure regularly?	50
<b>E</b>	<b>Technical information</b>	
	1. Error messages	52
	2. Customer service	53
	3. Technical data	53
	4. Original spare parts and accessories	54
	5. Metrological inspection	55
	6. Explanation of symbols	56
	7. Disposal	57

## Table of contents

<b>F</b>	<b>Maintaining the device</b>	58
<b>G</b>	<b>Warranty</b>	59

Thank you for choosing the visomat comfort XXL upper arm blood pressure monitor (hereinafter also referred to as the device).

These instructions are intended to help the user to use the unit safely and efficiently, and must be kept with the product and forwarded, if applicable. The unit must be used in accordance with the procedures contained in these instructions for use and must not be used for any other purposes.

It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 41.

This device is intended for measuring the blood pressure and pulse at the upper arm of persons over 12 years of age. It is designed for self-testing or for use by health care professionals.

The device has been clinically validated for patients with hypertension, hypotension, diabetes, pregnancy, pre-eclampsia, atherosclerosis, end-stage kidney disease and obesity, as well as for elderly patients.

Indications: Hypertension, (orthostatic), hypotension, white coat effect, masked hypertension, sleep apnoea syndrome, diabetes, organ damage.

Contraindications: Stents placed in the arteries of the arm to be measured, upper arm pain/bruising due to the frequency of high cuff pressure during measurement in patients with high blood pressure.

## 1. Important patient information

- The cuff must only be used on the upper arm and must never be used on other limbs.
- Taking blood pressure measurements on children requires specialist knowledge. Consult your doctor if you want to measure a child's blood pressure. Under no circumstances must the device be used on an infant.
- This blood pressure monitor has been specifically tested for use during pregnancy and by patients with pre-eclampsia. If you obtain abnormally high readings during pregnancy, repeat the measurement after approximately one hour. If the readings are still too high, consult your doctor or gynaecologist.
- Do not under any circumstances place the cuff over a critical area, e.g. a wound, aneurysm, etc. or on an arm with an arteriovenous shunt. Risk of injury! Any supply via an intravascular access point (infusion) or other medical monitoring devices could possibly be interrupted.
- Do not use the device without first consulting your doctor if you are having dialysis treatment or are taking anticoagulants, platelet aggregation inhibitors or steroids. These circumstances may cause internal bleeding.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other in terms of their mode of operation.
- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.

- The device contains small parts that could be swallowed by children. The air hose poses a strangulation hazard. Do not leave the device unattended with children.
- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axilla lymph nodes.
- The readings measured using this device do not provide a diagnosis. They are no substitute for visiting a doctor. Do not rely solely on the result obtained; always take the patient's other symptoms into account too. In case of doubt, please call a doctor or paramedics to help.
- Discuss the measured values with your doctor. Do not assess the measurement results yourself. **Never change the dose of medicines prescribed by your doctor.**
- Before carrying out your self-measurement, please pay attention to the chapter "Important usage information" on page 41.

## 2. Important technical information

- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight. See the storage and operating conditions in "Technical data" on page 53.
- Do not knock or drop the device.
- The device is not waterproof. Never immerse it in water or any other liquids. If liquid enters the device, this may result in serious damage and may cause malfunctions.
- The unit must only be operated with the cuff and accessories provided (see "Original spare parts and accessories" on page 54). Using different cuffs or accessories can lead to incorrect measure-

ment results. The warranty will be invalidated if the unit is damaged by third-party accessories!

- A consistently good power supply to the device is necessary for fault-free blood pressure measurements.
  - Only use long-life alkaline batteries (4 x 1.5 V, size AA/LR6).
  - Always replace all batteries simultaneously during battery replacement. Do not mix new and old batteries or batteries of different types. Do not use batteries beyond the specified expiry date.
  - When using the device with a mains adapter, please only use the visomat mains adapter, which has been specially tested for medical devices.
- Never simultaneously touch the device's live electrical parts and the patient.
- If the device is not to be operated for an extended period of time, please remove the batteries. In principle, all batteries can leak.
- Do not remove batteries from the device or separate it from the mains adapter when the device is switched on.
- Using this device in the vicinity of mobile phones, microwaves or other devices with strong electromagnetic fields can cause malfunctions and inaccurate measurements. During use, maintain a distance of at least 3.30 m from such devices.
- Never open or modify the device or the cuff (except for replacing the batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See "Technical data" on page 53.
- If the device has been stored at temperatures below 10°C or above 40°C, allow it to stand at room temperature for at least 2 hours before using it.

## Safety instructions

- The inflating and measuring procedure can be interrupted by pressing the Start/Stop button or by removing the cuff. The device then stops the inflation procedure and deflates the cuff.
- Before use, check that the device is operating reliably and is in proper condition. Do not use the device if it is damaged, as this can lead to injuries or inaccurate measurements.
- When using the power adapter, always make sure that the adapter and cable are not damaged.
- If the device malfunctions or is faulty, please refer to the troubleshooting guide starting on page 52, or contact customer service (see page 53).


## Operating the device

### 1. Device description

- |                                       |  |
|---------------------------------------|--|
| 1 Cuff connection                     | 6 Display                              |
| 2 Cuff                                | 7 Battery compartment                  |
| 3 Marking for arm circumference       | 8 Start/Stop button / user selection   |
| 4 Artery marking                      | 9 Memory button (memory recall button) |
| 5 Connecting socket for mains adapter |  |

## 2. Display

- |  |   |
|--|---|
| <b>1</b> Memory space or date/time                                     | <b>5</b> User                                     |
| <b>2</b> SYS = systole   | <b>6</b> Battery indicator                        |
| <b>3</b> DIA = diastole  | <b>7</b> Pulse signal display or irregular pulses |
| <b>4</b> PUL 1/min = pulse<br>Calculated pulse frequency<br>per minute | <b>8</b> WHO classification (page 47)             |
|  | <b>9</b> Memory code                              |

## 3. Important usage information

The measurement results of automated blood pressure monitors can be influenced by the measuring location, posture, previous exertion and general physical condition. Please pay attention to the usage information in order to obtain accurate readings.

- Do not consume alcohol, caffeine or smoke for at least one hour before measuring.
- Rest for at least 5 minutes before taking the measurement. Depending on the previous degree of exertion, you may even need to rest for up to an hour.
- Expose the upper arm. Clothing must not hinder the blood flow to or from the arm, as this will affect the blood pressure at the measuring point and may make the readings inaccurate. To avoid pinching, do not roll up the sleeves on shirts or blouses when taking a measurement. As long as the material lies flat against the arm, the measurement will not be affected.
- Do not move or talk during measuring.
- Breath slowly and deeply. Do not hold your breath.
- Check the display for irregular pulse signals. If required, repeat the measurement when conditions have improved.
- In order to exclude side differences and to obtain comparable measurement results, it is important to always take the measurement on the same arm. Ask your doctor which side it is advisable for you to measure.
- Always take your daily blood pressure measurements at around the same time.
- Blood pressure is not a fixed value. It can go up or down by more than 20 mmHg in patients within a few minutes.

#### 4. Initial operation of the device

Insert the batteries provided into the device.

If the device is to be operated using mains power, the cable connector of the mains adapter (not supplied) must be plugged into the socket on the rear of the device. The batteries are automatically switched off.

Please only use the visomat mains adapter. See “Original spare parts and accessories” on page 54.

#### 5. Inserting/replacing batteries, battery icon

- Remove the battery compartment cover on the underside of the device.
- Remove the old batteries from the device and insert the new ones. Please ensure correct polarity (marking in battery compartment).
- Close the battery compartment by clipping the battery cover back into the device.

When the “Empty battery 

If the batteries are completely flat, the empty battery symbol 

Once the batteries have been replaced, the readings will still be stored in the memory, but the date and time will need to be reset. To do this, follow the steps below.

#### 6. Activating and setting the date/time

The device has a date/time function. This function is activated on your device by default. If you do not need this function, you can disable it manually. This function can be reactivated later if required.

When the date/time function is activated, the measurement results are saved along with their date and time. This can make it easier to view the readings at a later point.

##### Activating and deactivating the date/time function

To activate the date/time function press and hold the START/STOP button. Then immediately press the MEMORY button while the full display is still visible. The function can be deactivated at any time using the same procedure.

When the date/time function is activated, the device continues to display the time even when it is switched off. The power consumption for this is extremely low.

##### Setting the date/time

Once the batteries have been inserted, the four-digit year format will appear on the display (Figure 1). Press the MEMORY button to set the year. Confirm the entry by pressing the START/STOP button.


Figure 1

You are now able to set the month. (Figure 2). Press the MEMORY button to set the month. Confirm the entry by pressing the START/STOP button.


Figure 2

Figure 3

Repeat these steps to enter the day (Figure 2) and the hours and the minutes (Figure 3).

After replacing the batteries, you will have to reset the date and time.

## Changing the date/time

To change the date and time, take the batteries out and wait until the display goes off. Then you can reset the date and time.

## 7. Attaching the cuff

Before applying the cuff, please check that the circumference of your upper arm is within the range specified on the cuff. An incorrect cuff size could result in inaccurate readings.

- Expose your upper arm.
- Push the cuff onto the arm until the lower edge of the cuff is 2 cm above the crook of the arm (Figure 1, Figure 2).
- When placed on the left arm, the air hose runs down the middle of the crook of the arm to the device. This automatically positions the 4 cm long artery marking centrally above the pulse position (Figure 3).
- When placed on the right arm, the cuff must be rotated to the left until the artery marking is on the pulse position. The air hose then runs along the inner side of the upper arm (Figure 4).
- Fasten the cuff, but make sure that it is not too tight. Tip:
  - Angle the arm slightly
  - Slightly tense the upper arm muscle
  - This increases the circumference of the arm slightly


Figure 1


Figure 2


Figure 3

Figure 4

Figure 5

- Now pull the free end of the cuff tight and close the Velcro fastener (Figure 5). If the cuff is looser at the lower end of your upper arm, fasten the cuff diagonally towards the shoulder to ensure that the cuff sits evenly across the length of the upper arm (Figure 6).

Figure 6

- Check that the metal bracket on the cuff is within the „Marking for arm circumference“ guide on the edge of the cuff (Figure 7).



Figure 7

- Connect the plug of the cuff to the cuff connection on the left-hand side of the device. Ensure that the plug is fully inserted into the device.

## 8. Body position during the measurement

- Sit at a table (preferably at the height of a dining table, not a coffee table).
- Sit with your back against the backrest of the chair.
- Rest your entire lower arm on the table with the palm facing upwards.
- Put your feet on the floor and do not cross your legs.

Ensure that the measuring point on the upper arm is approximately level with your heart. If the measuring point is lower than the heart, the measured values may be higher. If the measuring point is higher than the heart, the measured values may be lower.

## 9. Selecting the user

The device is designed for two users. The values for the two users are stored separately and can also be evaluated separately.

You may also use the two user memories to evaluate measurements taken in morning and those taken in the evening separately. To do this, carry out the measurement in the morning using user 1 and in the evening using user 2.

- Before taking the measurement, check that the correct user is displayed.
- To change the user, press and hold the START/STOP button until the display for the desired user number flashes on the screen. Press the MEMORY button to switch between User 1 and User 2. Confirm the entry by pressing the START/STOP button.
- Now carry out the measurement.

## 10. Measuring blood pressure

- Switch on the device by pressing the START/STOP button. All segments and icons will appear on the display (Figure 1).
- After five short beeps, the cuff will start to inflate automatically. The cuff will be inflated to the pressure required for measurement (Figure 2).
- The actual measurement process then begins as the cuff deflates. The ♥ symbol flashes in the display (Figure 3). In addition, beeps denote the pulse.

Figure 1

Figure 2

Figure 3

- A long beep denotes the end of measurement. The cuff is automatically deflated.
- The determined systole, diastole and pulse values are shown in the display (Figure 4).

Figure 4

The device switches itself off after approximately 3 minutes. You can also switch the device off before this by pressing the Start/Stop button.

## 11. Classification of the measured values (WHO)

The unit categorises the blood pressure readings according to the guidelines of the World Health Organisation (WHO). You can use the coloured fields to read the categorisation of each reading:

- Upper dot (Red) signifies high pressure\*
- Middle dot (Yellow) signifies borderline value\*
- Lower dot (Green) signifies normal value

WHO 2003	Systolic pressure = Upper value mmHg**	Diastolic pressure = Lower value mmHg**
Red	from 140*	from 90*
Yellow	120 to 139	80 to 89
Green	below 120	below 80

\* It is sufficient for one of the values to be raised.  
\*\* millimetres of mercury

The WHO classification is stored along with the measured values and both can be recalled from the memory.

Depending on age, weight and general health, the blood pressure values can be different. Only a doctor can determine the right blood pressure range for you and assess whether your blood pressure has

reached a dangerous level for you. Discuss your blood pressure values with your doctor.

**Never change the dose of medicines prescribed by your doctor!**

## 12. Irregular pulses and cardiac arrhythmia

If the device detects irregular pulses during the measurement, the  symbol is displayed once the measurement is complete.

This can be triggered by cardiac arrhythmias, movement, talking or even by deep breathing. The symbol is stored with the respective measurement.

A more frequent appearance of the symbol may indicate cardiac arrhythmia and must be discussed with your doctor! Measurement results accompanied by irregular pulses should be viewed as critical and repeated under more favourable conditions.

The device must not be used in place of a cardiac exam; it is only intended for early detection of an irregular pulse.

## 13. Using the memory

The device stores the measured values in two memories (one per user). Each memory can store up to 60 readings and the average value (A).

The measured results are automatically saved in the memory. If more than 60 measured values are saved, the oldest measured value (no. 60) is deleted so that the latest value (no. 1) can be recorded.

### Data retrieval

Set the user whose data you wish to retrieve (see page 46).

To retrieve data, press the MEMORY button. The average of the saved results will be indicated by an "A", followed by the number of saved readings from which the average was calculated (Figure 1).

Figure 1

Press the Memory button again to display the most recently measured values (Figure 2). If the date/time function is activated, the memory number (Figure 2), the date (Figure 3) and the time (Figure 4) are displayed alternately.

Figure 2

Figure 3

Figure 4

Press the Memory button repeatedly to recall further measured values.

### Deleting data

To delete individual measured values, recall the measured value you wish to delete.

Press and hold the MEMORY button again until the reading flashes several times and then goes out.


Figure 5

To delete the entire memory, retrieve the respective average value (A), press the Memory button again and hold down until the average value firstly flashes and is then replaced with dashes (Figure 5).

## **D** What you should know about blood pressure

### **1. The systolic and diastolic blood pressure value**

The circulatory system has the important task of supplying all organs and tissues in the body with sufficient blood, and of removing metabolites. The heart does this by contracting and expanding at a regular rhythm of around 60-80 times per minute. The pressure created by the blood flow on the artery walls when your heart contracts is referred to as systole. The pressure in the subsequent relaxation phase when the heart fills with blood again is referred to as diastole. In your daily measurement, you will measure both values.

### **2. Why you measure different values**

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant.

Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure.

This explains why values measured at the doctor are often higher than those you obtain at home in your usual environment.

### **3. Why should you measure blood pressure regularly?**

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up.

## **D** What you should know about blood pressure

One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is only possible when measurements are taken regularly. Discuss the measured values with your doctor.

## 1. Error messages

Error	Possible cause	Remedy
Display Err 1	The systolic pressure was measured, but then the cuff pressure dropped to less than 20 mmHg. This situation occurs if the tube was disconnected after the systolic blood pressure was measured. Other possible cause: No pulse was detected.	Severe leakage while the measurement is being taken. Check the connection. If the connection is OK, please send the device to us.
Display Err 2	Unnatural pressure pulses will affect the result. Cause: The arm moved while the measurement was being taken (artefact).	Keep still.
Display Err 3	The cuff takes a long time to inflate. The cuff is not sitting correctly or the connection to the tube is not airtight.	If this is not the case, please send the device to us. If there are multiple issues, please call customer services.
Display Err 5	The device has detected an implausible difference between the systolic and the diastolic pressure. Check that the readings are plausible. The device will question readings such as 160/140.	Take another measurement and follow the instructions carefully. If you continue to obtain abnormal readings, contact your doctor.
Display HI	The inflation pressure exceeds 300 mmHg or the pulse is above 200. The device can only display a pulse of 40–200. Possible arrhythmia, respiratory arrhythmia.	Arrhythmia — rest for 5 minutes, then repeat the measurement. Respiratory arrhythmia — rest for 5 minutes, then repeat the measurement. Take shallow breaths. If necessary, check with your doctor.
Display LO	Pulse reading is below 40. The device can only display a pulse of 40–200.	Repeat the measurement after resting for 3–5 minutes. Intermittent pulse sounds? Listen carefully to the beep and check with your doctor if necessary.

Error	Possible cause	Remedy
	Irregular pulse (cardiac arrhythmia, movement during the measurement, tremors, shaking, deep breathing, etc.)	Repeat the measurement after resting for 3–5 minutes. Always discuss anything that could have affected the result (depending on the severity of the arrhythmia) with your doctor.
Readings are too high	Did the patient rest for the specified amount of time before the measurement was taken?	Rest for approximately 3–5 minutes, then repeat the measurement. - Do not move your arm - Do not speak
Abnormal readings	Movement or talking while the measurement is being taken. Patient did not rest beforehand.	Check the conditions (see "Important usage instructions" on page 41) and repeat the measurement.

## 2. Customer service

Device repairs may only be performed by the manufacturer or an expressly authorised body. Please contact:

UEBE Medical GmbH info@uebe.com  
 Bürgermeister-Kuhn-Straße 22 www.uebe.com  
 97900 Kulsheim, Germany

## 3. Technical data

Model:	visomat comfort XXL
Size:	L = 162 mm x W = 110 mm x H = 80 mm ± 2 mm
Weight:	approx. 310 g without batteries
Display:	LCD display (liquid crystal display) for measured values and check displays
Memory:	2 memories each capable of storing 60 readings
Measurement method:	Oscillometric determination of systole, diastole and pulse
Reference method of the clinical test:	Auscultatory measurement

Pressure display range:	0 - 299 mmHg
Measurement range:	Blood pressure measurement: 30 - 280 mmHg Pulse measurement: 40 - 200 pulses/minute
Measurement precision:	Blood pressure measurement (cuff pressure): $\pm 3$ mmHg, Pulse rate: $\pm 5\%$
Power supply:	4 x 1.5 V AA LR6 alkaline manganese batteries, durability: at least 400 measurements Optional: visomat mains adapter, output 6 V DC, 500 mA
Power consumption:	max. 4 W
Cuff:	Cuff type WXL for an arm circumference of 32 - 52 cm, applied part type BF
Operating conditions:	Ambient temperature 10 °C to 40 °C, relative humidity 15 % to 85 % non-condensing, air pressure 700 to 1060 hPa
Conditions for storage and transport:	Ambient temperature -20 °C to +50 °C Relative humidity 15 % to 85 % $\pm 2\%$ , non-condensing
Time to adjust from minimum/maximum storage temperature to operating temperature:	2 hours
Inflation and deflation:	Regulated automatically
IP rating:	IP 20: Protected against solid foreign bodies with a diameter from 12.5 mm, no protection against water.
Protection against electric shock:	Internal power supply, applied part type BF (cuff)
Expected service life:	Device: 5 years or 10,000 measurements Accessories: 2 years
Operating mode:	Continuous operation
Classification:	Internal power supply using batteries.

#### 4. Original spare parts and accessories

The following original spare parts or accessories are available through specialist retailers:

EN-54

- Universal cuff 32 - 52 cm, type WXL  
Part no. 2407601, PZN-16086469
- Cuff 32 - 42 cm, type L2  
Part no. 2402502, PZN-07333408
- Cuff 22 - 32 cm, type M2  
Part no. 2402501, PZN-03088490
- visomat mains adapter  
Part no. 2400020, PZN-11870684

Subject to technical modifications.

#### 5. Metrological inspection

Generally, it is recommended that a metrological inspection is performed every two years. However, professional users in Germany are obligated to do so in accordance with the "Medical Products Operator Ordinance" (Medizinprodukte-Betreiberverordnung).

This can be performed either by UEBE Medical GmbH, an authority responsible for metrology or an authorised maintenance service. Please refer to your national regulations.

Upon request, responsible authorities or authorised maintenance services receive a „Test instruction for metrological inspection“ from the manufacturer.

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

Important: No modifications, e.g. opening the device (except to replace the batteries), may be made to this device without the manufacturer's permission.

EN-55

## 6. Explanation of symbols

 **0123** This product complies with Council Directive 93/42/EEC of 5 September 2007 concerning medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock TYPE BF

Consult the instructions for use

Keep dry

Device serial number

Reference number = item number

**IP20**

Protected against solid foreign bodies with a diameter from 12.5 mm, no protection against water.

Manufacturer

Polarity inside plus

Cuff connection

### Mains adapter information

When using the optionally available mains adapter:

Protection class II (double insulation)

For indoor use only

Direct current

## 7. Disposal

Technical equipment and batteries do not belong in household waste. They must be disposed of at appropriate collection and disposal points.

## Maintaining the device

- Keep the device clean. Check for cleanliness after use. Please use a soft, dry cloth for cleaning. Do not use benzene, thinners or other strong solvents.
- The cuff may absorb sweat and other fluids. Inspect the cuff for marks or discolourations after each use. Clean the device using soapy water no warmer than 30 °C and carefully wipe the surface. Do not scrub or machine-wash. Air-dry thoroughly.
- Make sure that liquid cannot get into the air hose.
- Do not apply force when folding or bending the cuff.
- To separate the cuff from the main device, please do not pull the air hose, but grip the air connector and gently pull it off.

## Warranty

The device has been manufactured and tested with all due care. Nevertheless, in the event of defects upon delivery, we provide a warranty under the following conditions:

During the warranty period of 5 years from the date of purchase, we remedy such defects at our discretion and at our expense in our plant through repair work or replacement delivery of a defect-free device.

The warranty does not cover normal wear and tear of wearing parts or damage caused by failure to observe the instructions for use, improper handling (e.g. unsuitable power sources, breakage, leaking batteries) and/or disassembly of the device by the buyer. Furthermore, the warranty does not constitute grounds for asserting claims for damages against us.

Warranty claims can only be asserted during the warranty period and upon presentation of the proof of purchase. In the case of a warranty claim, the device along with the proof of purchase and a description of the complaint must be sent to:

UEBE Medical GmbH  
Service-Center  
Zum Schlag 18  
97877 Wertheim, Germany

The cost of returning the device to our factory shall be borne by the sender. Complaints that are sent back without prepayment will not be accepted by UEBE.

The statutory claims and rights of the buyer against the seller (for example, defect claims, producer liability) are not restricted by this warranty.

**Please note: In the case of a warranty claim, please make sure to enclose the proof of purchase.**

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© Copyright 2021 UEBE Medical GmbH